



29 Palms Laboratory  
 47-250 Dillon Road  
 Coachella, Ca 92236  
 Phone: 760-398-6767  
 Fax: 760-398-0046

Title: Portable Turbidimeter  
 Number: PP0010  
 Release Date: April 16, 2003  
 Revision Date: April 16, 2003  
 Version: 1.0

**DOCUMENT TYPE:** Standard Operating Procedure

**DOCUMENT CLASS:** Physical Property Procedure

**TITLE:** Portable Turbidimeter, Model 2100  
 US EPA, Method 180.1

**PREPARED BY:** Alison Millar, Environmental Technician II \_\_\_\_\_

**REVISED BY:** Alison Millar, Environmental Technician II \_\_\_\_\_

**REVIEWED BY:** Stephanie Hull, Environmental Scientist II \_\_\_\_\_

**APPROVED BY:** Marshall K. Cheung, Ph.D., Laboratory Director \_\_\_\_\_

## Table of Contents

1.	Scope and Application.....	3
2.	Summary of Method.....	3
3.	Comments.....	3
4.	Sample Handling and Preservation .....	3
5.	Interferences .....	3
6.	Apparatus.....	3
7.	Reagents .....	4
8.	Procedure.....	4
9.	Calculation.....	5
10.	Results .....	5
11.	Bibliography.....	5
12.	Calibration Curve .....	6



29 Palms Laboratory  
47-250 Dillon Road  
Coachella, Ca 92236  
Phone: 760-398-6767  
Fax: 760-398-0046

Title: Portable Turbidimeter  
Number: PP0010  
Release Date: April 16, 2003  
Revision Date: April 16, 2003  
Version: 1.0

<b><i>Document No.:</i></b>	<b><i>PP0010-005</i></b>
<b><i>Copy provided to:</i></b>	<b>Marshall K. Cheung</b>
<b><i>Title:</i></b>	Laboratory Director 29 Palms Laboratory 47-250 Dillon Road Coachella, CA 92236
<b><i>Copy provided by:</i></b>	<b>Alison Millar</b>
<b><i>Title:</i></b>	Environmental Technician II
<b><i>Date:</i></b>	July 28, 2004



## 1. Scope and Application

- 1.1. This method is applicable to the analysis of drinking, surface and saline waters in the range of 0 to 1000 nephelometric turbidity units (NTU).
- 1.2. More turbid samples should be diluted prior to analysis.

## 2. Summary of Method

- 2.1. The method is based upon a comparison of the intensity of light scattered by the sample under defined conditions with the intensity of light scattered by a standard reference suspension under the same conditions.
- 2.2. The higher the intensity of scattered light, the higher the turbidity.
- 2.3. StablCal standard suspensions are used to create a calibration curve.
- 2.4. Samples are transferred to sample cells that are carefully wiped free of dust and moisture.
- 2.5. Cuvettes are inserted into the chamber with a consistent orientation.
- 2.6. Turbidity units are recorded according to results. (See Results, p.5)

## 3. Comments

- 3.1. NTU's are considered comparable to the previously reported Formazin Turbidity Units (FTU) and Jackson Turbidity Units (JTU).
- 3.2. Dark glass sample cell should be handled by the securely fastened lid.
- 3.3. Store StablCal standards at room temperature, preferably in drawer labeled Turbidity in 29 Palms Laboratory.
- 3.4. Keep lit closed when not in use.
- 3.5. Make sure that there are no air bubbles in sample cell before placing into Turbidimeter.
- 3.6. Wipe sample cell well with Kimwipe to remove any fingerprints, and apply a drop of silicone oil to exterior of cell to mask any minor imperfections in the surface.
- 3.7. A Formazin calibration should be done once every 3 months

## 4. Sample Handling and Preservation

- 4.1. Preservation of the sample is not practical; analysis should begin as soon as possible.
- 4.2. Refrigeration or icing to 4°C, to minimize microbiological decomposition of solids is recommended.

## 5. Interferences

- 5.1. Coarse floating debris that settles out rapidly will give unstable low readings.
- 5.2. Presence of color in sample water due to dissolved substances that absorb light will cause turbidities to be low.
- 5.3. Air bubbles will cause higher turbidity readings.

## 6. Apparatus

- 6.1. Portable Turbidimeter, Model 2100P
- 6.2. Silicone Oil
- 6.3. Oiling cloth
- 6.4. Quick reference card
- 6.5. Four AA batteries



## 7. Reagents

- 7.1. Reagent Water (NanoPure)
- 7.2. StablCal Turbidity Standards:  $\leq 0.1$ -NTU, 20-NTU, 100-NTU, 800-NTU

## 8. Procedure

- 8.1. Turn instrument on
- 8.2. Calibrate the meter.
  - 8.2.1. Insert the “0.1” turbidity sample cell calibration standard into the cell compartment by aligning the orientation mark on the cell with the mark on the front of the cell compartment.
  - 8.2.2. Close shield.
  - 8.2.3. Press CAL and the CAL and S 0 icons will be displayed with the “0” will be flashing.
  - 8.2.4. Press READ on keypad.
    - 8.2.4.1. The instrument will count from 60 to 0, read the blank and use it to calculate a correction factor for the next standard measurement. The display will automatically increment to the next standard.
  - 8.2.5. Remove sample cell.
  - 8.2.6. The display will show S 1 with the “1” flashing and the 20 NTU value.
  - 8.2.7. Insert the “20” cell into the well compartment and close cover.
  - 8.2.8. Press READ.
    - 8.2.8.1. The instrument will count from 60 to 0, read the blank and use it to calculate a correction factor for the next standard measurement. The display will automatically increment to the next standard.
    - 8.2.8.2. If the value is incorrect, edit the value by pressing the → key to scroll to the correct number.
    - 8.2.8.3. Repeat sample reading with new 20 NTU standard.
  - 8.2.9. Remove sample cell.
    - 8.2.9.1. The display will show the S 2 with the “2” flashing.
  - 8.2.10. Insert “100” NTU standard sample into the well compartment.
  - 8.2.11. Press READ.
    - 8.2.11.1. The instrument will count from 60 to 0, read the blank and use it to calculate a correction factor for the next standard measurement. The display will automatically increment to the next standard.
    - 8.2.11.2. If the value is incorrect, edit the value by pressing the → key to scroll to the correct number.
  - 8.2.12. Remove calibration sample cell.
    - 8.2.12.1. The display will show the S 3 with the “3” flashing and 800 NTU for the next sample.
  - 8.2.13. Place 800 NTU calibration sample cell into well compartment.
  - 8.2.14. Press READ.
    - 8.2.14.1. The instrument will count from 60 to 0.
    - 8.2.14.2. The display will then increment back to SO.
  - 8.2.15. Remove calibration sample cell.
  - 8.2.16. Press CAL to accept the calibration.
- 8.3. Perform Calibration Curve
  - 8.3.1. Place a known standard in cell compartment



**29 Palms Laboratory**  
47-250 Dillon Road  
Coachella, Ca 92236  
Phone: 760-398-6767  
Fax: 760-398-0046

**Title: Portable Turbidimeter**  
**Number: PP0010**  
**Release Date: April 16, 2003**  
**Revision Date: April 16, 2003**  
**Version: 1.0**

- 8.3.2. Press READ
- 8.3.3. Record Value
- 8.3.4. Place different known standard into cell well.
- 8.3.5. Press READ
- 8.3.6. Record Value
- 8.3.7. Repeat 8.3.1 three more times for a total of five recorded values.
  - 8.3.7.1. Proceed to section 9.0 for calculations
- 8.4. Take sample readings.
  - 8.4.1. Fill sample cell to white horizontal line with sample.
  - 8.4.2. Wipe sample cell well with Kimwipe and lightly oil with silicone and wipe off.
  - 8.4.3. Insert sample cell with vertical line forward.
  - 8.4.4. Close lid.
  - 8.4.5. Press READ
  - 8.4.6. Record value of sample in NTU's.

## 9. Calculation

- 9.1. Using Excel graph the difference between the expected NTU value and the recorded actual value to find the slope. (See section 12.0 for Calibration Curve)
- 9.2. If created diluted sample.
  - 9.2.1. Multiply sample readings by appropriate dilution to obtain final reading.
  - 9.2.2. Sample calculation: If sample is diluted in half (1:1), the turbidity result would be the value displayed multiplied by two.
    - 9.2.2.1. Diluted sample turbidity reading: 5.4.
    - 9.2.2.2. Turbidity (NTU) =  $5.4 \times 2 = 10.8$  NTU
    - 9.2.2.3. Reported as 11 NTU (See Results p.5)

## 10. Results

- 10.1. Report results as follows:


<u>NTU</u>	<u>Record to nearest:</u>
0.0-1.0	0.05
1-10	0.1
10-40	1
40-100	5
100-400	10
400-1000	50

## 11. Bibliography

- 11.1. Hach Company Manual for Portable Turbidimeter, Model 2100P (Cat. No. 46500-88)
- 11.2. EPA Method 180.1
- 11.3. Hach Quick Reference Card for Portable Turbidimeter.



## 12. Calibration Curve



Turbidity Expected (NTU)	Measured Turbidity (NTU)	RF	
0.00	0.00	-	Average RF =0.964 Std. Dev. =0.0492 Relative % Std. Dev. =5.10  R <sup>2</sup> =0.997 m =0.913 y-intercept =0.00
2.00	2.00	1.000	
5.00	4.97	0.994	
10.00	9.68	0.968	
20.00	17.86	0.893	